

DEXXUM T BEYOND DEXXUM T

Quantum

New fan beam DXA



- Central DXA(Dual energy X-ray Absorptiometry)
- Fan beam
- Scan site : AP spine, Femur(Dual femur), Forearm
- Scan area : 580×480 mm
- Scan time : AP spine - 36Sec.
Femur - 27 Sec.
Forearm - 24Sec.
- Automatic real one-scan
- Reproducibility : $\leq 1.0\%$ C.V.
- Measured parameter : BMD, BMC, BMI, T-score, Z-score, Area
- FRAX / Color mapping / Trend report / DICOM & PACS
- Dimension : $2000 \times 800 \times 1213$ mm(Standard)
 $1900 \times 800 \times 1213$ mm(Medium)
 $1850 \times 800 \times 1213$ mm(Compact)
- Table height : 650mm
- Weight : 130.5Kg
- Power consumption : 110VAC / 220VAC (+/- 10%)

DEXXUM T *Quantum*


Remote Easy Plug Easy Connect Easy Clean Easy Install Easy Use

[illegible]

Femur



spine

[illegible]

OsteoSys

9F, 903 JnK Digital Tower, 111 Digital-ro 26, Guro-gu, Seoul, KoreaTel. +82.2.6124.5900 Fax. +82.2.6124.5958
www.osteosys.com